InBody BWA

Connecting and Registering LB Web User Manual

CONNECT

Connect to internet by either Wi-Fi or LAN

WI-FI (WIRELESS)


Learn how to connect InBody to Wi-Fi

Connect InBody – Wi-Fi



Press **[Administrator Menu]**.


Connect InBody – Wi-Fi



Input your 4 digit password and press **[Enter]**.


Note: By default the access code will be **[0 0 0 0]**.

Connect InBody – Wi-Fi




Use the right arrow **[>]** then press **[16. Internet Options]**.

Connect InBody – Wi-Fi




Press **[O]**.

Connect InBody – Wi-Fi



Check the **[Wi-Fi (Wireless)]** box and press **[Next]**.

Connect InBody – Wi-Fi



Check the Wi-Fi connection box you wish to connect to and press **[Connect]**.

Note: If you do not see your wireless signal press the **[Refresh]** button or scroll down.

Connect InBody – Wi-Fi




Input the Wi-Fi password and **[Enter]**.

Connect InBody – Wi-Fi



Please wait...

Connect InBody – Wi-Fi



Once you have receive an IP Address, you are now connected.

LAN (Wired)


Learn how to connect your InBody via LAN

Connect InBody – LAN




Press **[Administrator Menu]** input the access code then **[Enter]**.

Connect InBody – LAN




Use the right arrow **[>]** then press **[16. Internet Options]**.

Connect InBody – LAN




Press **[O]**.

Connect InBody – LAN



Check **[LAN (Wired)]** box and press **[Next]**.

Connect InBody – LAN




Press **[Next]**.

Connect InBody – LAN



Please wait...

Connect InBody – LAN



Once you receive an IP Address, you are now connected.

CLOUD SERVICE


Learn how to connect your InBody to the Cloud Service

InBody – Cloud Service



Press **[Administrator Menu]** input the access code then **[Enter]**.

InBody – Cloud Service




Select **[04. Cloud Service]**.

InBody – Cloud Service



Please wait while the InBody checks for a stable internet connection.

InBody – Cloud Service



Press the top left corner of the LCD screen where **[04. Cloud Service]** is located five times.

Note: If this screen does not appear for you skip to page 30.


InBody – Cloud Service



Enter in **EXACTLY** as followed **https://usa.lookinbody.com/InBodyAPI/** and **[Enter]**.


Note: Every letter is case sensitive and please make sure to enter in the " / " at the end of the URL.

InBody – Cloud Service



Press **[Confirm]**.

InBody – Cloud Service




Select **[04. Cloud Service]** once more.

InBody – Cloud Service




Please wait while the InBody checks for a stable internet connection.

InBody – Cloud Service




Select the **O**.

InBody – Cloud Service



Check the **[Lookin'Body Website Account Login]** box.


InBody – Cloud Service



Enter your ID and Equip Password then select **Login**.


Note: If you do have not received any login information for your InBody, please visit our Download Center and fill out the Web application form.

InBody – Cloud Service




Select **Save** to confirm the changes.

InBody – Cloud Service



Press Exit to return to the main screen.

InBody – Cloud Service



You are now ready to use the Lookin'Body Web!